FORTGESCHRITTENEN PRAKTIKUM II

Optical pumping

11.04.2016

Benjamin Winkelmann Peter Spalthoff

Tutor: Dominik Schomas

Contents

Li	st of Figures
1	Goal of the experiment
2	Physical principles
	2.1 Hyperfine structure and Zeeman splitting
	2.2 Optical pumping
	2.3 Relaxation processes
	2.4 Larmor precession of the spin
	2.5 The laser diode
3	References
4	References

List of Figures

2.1	Hyperfine structure of Rubidium	4
2.2	Hyperfine structure energies	•
2.3	Optical pumping	4

1 Goal of the experiment

In this experiment, the process of optical pumping will be used to precisely measure properties of Rubidium atoms such as the hyperfine constant A via absorption measurements. In addition to that, relaxation times of the induced pumped states as well as external magnetic fields will also be measured by observing the effect of magnetic fields, applied through Helmholtz coils, high frequency radio waves and variations in the laser intensity.

2 Physical principles

2.1 Hyperfine structure and Zeeman splitting

This section is based on the detailed elaborations in [1].

The fine structure levels of the atomic spectrum, which splits the basic levels into sublevels due to spin-orbit interaction, can be shown to be split into even finer levels, whose energetic distances are roughly three orders of magnitude smaller than those of the fine structure. This is called the 'hyperfine structure' and is mainly caused by the interaction of the nuclear magnetic dipole and quadrupole moment and the magnetic field of the shell electrons. Its structure for the two Rubidium isotopes that are used in this experiment can be seen in figure 2.1.

As the nucleus is charged and, expressed as the nuclear spin \vec{I} , has angular momentum, it also has a magnetic moment, which is $\vec{\mu}_I = \frac{g_I \mu_K}{\hbar} \vec{I}$, where g_I is the g-factor of the nucleus and μ_K is the nuclear magneton.

With the total angular momentum of the electrons \vec{J} , the total angular momentum of the atom can be written as

$$\vec{F} = \vec{J} + \vec{I}, \qquad |I - J| \le F \le I + J \tag{2.1}$$

The energy difference between hyperfine structure levels can then shown to be

$$\Delta E_{HFS} = -\vec{\mu}_I \cdot \vec{B}_J = \frac{A}{2} (F(F+1) - J(F+1) - I(I+1))$$
 (2.2)

where $A = \frac{g_I \mu_K B_J}{\sqrt{J(J+1)}}$ is the hyperfine constant. Neighboring levels thus have an energy difference of

$$\Delta E_{HFS}(F+1) - \Delta E_{HFS}(F) = A(F+1) \tag{2.3}$$

This structure for the rubidium isotopes used in this experiment can be seen in figure 2.2.




Figure 2.1: The hyperfine structure of the two isotopes of Rubidium used in the experiment. The hyperfine levels in turn are split due to the Zeeman effect caused by an external field of $B=1~\mathrm{G}$. [1]




Figure 2.2: Energies of a set of hyperfine levels of the two isotopes used in the experiment. [1]

These hyperfine levels can in turn be split into 2(F+1) sub-levels in the presence of an external magnetic field. The according quantum number is $-F \leq m_F \leq F$. As long as the magnetic field is weaker than the spin-orbit coupling or, in other terms, $g_J \mu_B B_0 \ll A$, this is called the Zeeman effect. The effect for larger fields, where the spin-orbit coupling is disrupted, is called Paschen-Back effect.

For the Zeeman effect, the energy difference of the levels is

$$\Delta E_{Zeeman} = \frac{g_J}{2(I + \frac{1}{2})} \mu_B B_0 \tag{2.4}$$

2.2 Optical pumping

In general, pumping refers to constantly transferring electrons into higher energy levels until significantly more electrons are in the higher than in the lower state. This is called population inversion.

In the case of this experiment, this is done using a laser diode. As the goal is to examine magnetic fields using the Zeeman splitting, a way must be found to create population inversion within a single non-degenerate hyperfine structure level. Normally, electrons are equally distributed between said levels.

The selection rules for transitions

$$\Delta F = 0, \pm 1 \qquad (F = 0 \leftrightarrow F = 0)$$

$$\Delta m_F = 0, \pm 1 \qquad (2.5)$$

allow for a convenient way to change that. If only σ^+ -polarized light is used, only transitions with $\Delta m_F = +1$ are caused. Since the following decay is random within the bounds of the transition rules, the laser will pump all electrons into the $^2S_{1/2}$ state with

 $m_F = +2$, F = 2 for ⁸⁷Rb and $m_F = +3$, F = 3 for ⁸⁵Rb. Figure 2.3 illustrates this for two exemplary transitions.

Mathematically, this process can be described as

$$\left(\frac{dn}{dt}\right)_{P} = \frac{N-n}{T_{P}} \tag{2.6}$$

where n is difference of the levels in the two-level system, N the overall number of atoms in the system and T_P the characteristic pumping time of the system.




Figure 2.3: Optical pumping for ⁸⁷Rb. The σ^+ polarized light can only cause transitions with $\Delta m_F = +1$, thus achieving the desired pumping effect.

2.3 Relaxation processes

The desired pumping effect is counteracted mainly by three relaxation effects:

Diffusion to the wall: Upon hitting the glass containment, the rubidium atoms may lose their polarization. This process is inhibited by a buffer gas, which limits the mean free path of the atoms.

Collisions with the buffer gas: The rubidium atoms may use their polarization due to collisions with the buffer gas. The cross section of this event depend highly on what kind of gas is used. Best results are achieved with noble gases - in the case of this experiment, krypton was used.

Spin exchange When rubidium atoms collide, they may interchange their spins. While the overall polarization is preserved, the decoupling of nuclear and electron spins lead to a faster relaxation time. Much more detailed elaborations can be found in [2].

Overall, the relaxation can be described by the following differential equation

$$\left(\frac{dn}{dt}\right)_{R} = -\frac{n}{T_{R}}\tag{2.7}$$

where T_R is the characteristic relaxation time. A value of $T_R^{theo} = [6.5]ms$ is given in [1]. The overall process of polarization orientation is thus the sum of equations 2.6 and 2.7:

$$\left(\frac{dn}{dt}\right)_{O} = \left(\frac{dn}{dt}\right)_{P} + \left(\frac{dn}{dt}\right)_{R} = \frac{N}{T_{P}} - n\left(\frac{1}{T_{P}} + \frac{1}{T_{R}}\right)$$
(2.8)

The solution of this equation is an exponential

$$n(t) \propto e^{-\frac{t}{\tau}} \tag{2.9}$$

where $\tau = \frac{1}{T_P} + \frac{1}{T_R}$.

2.4 Larmor precession of the spin

If the ensemble is polarized along a certain magnetic field and one component of said field is suddenly set to zero, the polarization precesses around the remaining field. The precession frequency is

$$f_L = \frac{g_F \mu_B}{b} \cdot B =: \alpha \cdot B \tag{2.10}$$

where g_F are the Landé factors for the rubidium isotopes, quantified by Baur [1] as $g_F(^{85}Rb) = 1/3$ and $g_F(^{87}Rb) = 1/2$. The proportionality constant between the frequency and the remaining magnetic field thus is

$$\alpha(^{85}Rb) = 4.665 \,^{\text{kHz}}/\mu\text{T} \quad \alpha(^{87}Rb) = 6.998 \,^{\text{kHz}} \,^{\text{T}}$$
 (2.11)

2.5 The laser diode

3 References

4 References

- [1] Baur, Clemens. "Einrichtung des Versuchs Optisches Pumpen mit Laserdioden." Zulassungsarbeit, Freiburg 1997
- [2] Happer, William. "Optical pumping." Reviews of Modern Physics 44.2 (Apr. 1972): 170-238